

For Research Use Only

CELSR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83207-2-PBS

Basic Information

Catalog Number: 83207-2-PBS	GenBank Accession Number: NM_014246	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9620	CloneNo.: 240110D7
Source: Rabbit	UNIPROT ID: Q9NYQ6	
Isotype: IgG	Full Name: cadherin, EGF LAG seven-pass G-type receptor 1 (flamingo homolog, Drosophila)	
	Calculated MW: 329 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00162-2, CELSR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83207-2-PBS. Detection antibody: 83207-3-PBS. Standard: SY00794. Range: 0.156-10 ng/mL