

For Research Use Only

# CELSR1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83207-3-PBS

## Basic Information

Catalog Number:

83207-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_014246

GeneID (NCBI):

9620

UNIPROT ID:

Q9NYQ6

Full Name:

cadherin, EGF LAG seven-pass G-type receptor 1 (flamingo homolog, Drosophila)

Calculated MW:

329 kDa

Purification Method:

Protein A purification

CloneNo.:

240110F7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00162-1, CELSR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83207-3-PBS. Detection antibody: 83207-1-PBS. Standard: SY00794. Range: 0.156-10 ng/mL



Cytometric bead array standard curve of MP00162-2, CELSR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83207-2-PBS. Detection antibody: 83207-3-PBS. Standard: SY00794. Range: 0.156-10 ng/mL