

For Research Use Only

CEP350-Specific Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: **83208-2-PBS**

Basic Information

Catalog Number: 83208-2-PBS	GenBank Accession Number: NM_014810	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9857	CloneNo.: 230320D4
Source: Rabbit	UNIPROT ID: Q5VT06	
Isotype: IgG	Full Name: centrosomal protein 350kDa	
	Calculated MW: 351 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00178-1, CEP350-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83208-3-PBS. Detection antibody: 83208-2-PBS. Standard: SY00642. Range: 0.313-20 ng/mL



Cytometric bead array standard curve ofMP00178-2, CEP350-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83208-2-PBS. Detection antibody: 83208-1-PBS. Standard: SY00642. Range: 1.25-20 ng/mL