

For Research Use Only

IGDCC3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83211-2-PBS

Basic Information

Catalog Number:

83211-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10276

GenBank Accession Number:

BC067107

GeneID (NCBI):

9543

UNIPROT ID:

Q8IVU1

Full Name:

immunoglobulin superfamily, DCC subclass, member 3

Calculated MW:

814 aa, 87 kDa

Purification Method:

Protein A purification

CloneNo.:

230449D8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00192-1, IGDC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83211-2-PBS. Detection antibody: 83211-1-PBS. Standard: Ag10276. Range: 1.25-20 ng/mL.