

For Research Use Only

NDUFB8 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83216-2-PBS

Basic Information

| | | |
|--|---|---|
| Catalog Number: 83216-2-PBS | GenBank Accession Number: BC000466 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 4714 | CloneNo.: 230510F6 |
| Source: Rabbit | UNIPROT ID: O95169 | |
| Isotype: IgG | Full Name: NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 8, 19kDa | |
| Immunogen Catalog Number: AG6495 | Calculated MW: 22 kDa | |

Applications

Tested Applications:
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00204-1, Human NDUFB8 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83216-1-PBS. Detection antibody: 83216-2-PBS. Standard: Ag6495. Range: 2.5-20 ng/mL.