

For Research Use Only

NUDCD2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83217-4-PBS

Basic Information

Catalog Number: 83217-4-PBS	GenBank Accession Number: BC017934	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 134492	CloneNo.: 230513F10
Source: Rabbit	UNIPROT ID: Q8WVJ2	
Isotype: IgG	Full Name: NudC domain containing 2	
Immunogen Catalog Number: AG15579	Calculated MW: 157 aa, 18 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00205-2, NUDCD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83217-2-PBS. Detection antibody: 83217-4-PBS. Standard: Ag15579. Range: 2.5-40 ng/mL.



Bi-layer interferometry (BLI) kinetic assays of 83217-4-RR against Human NUDCD2 were performed. The affinity constant is 76.2 nM.