

For Research Use Only

Rat KIM-1/HAVCR1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83221-4-PBS

Basic Information

Catalog Number: 83221-4-PBS	GenBank Accession Number: NM_173149.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 286934	CloneNo.: 230528G9
Source: Rabbit	UNIPROT ID: O54947	
Isotype: IgG	Full Name: hepatitis A virus cellular receptor 1	
	Calculated MW: 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00213-2, Rat KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83221-2-PBS. Detection antibody: 83221-4-PBS. Standard: Eg0597. Range: 0.625-20 ng/mL.