

For Research Use Only

PLTP Recombinant antibody, PBS Only (Detector)



Catalog Number: 83222-2-PBS

Basic Information

Catalog Number: 83222-2-PBS	GenBank Accession Number: BC019898	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5360	CloneNo.: 230542F9
Source: Rabbit	UNIPROT ID: P55058	Affinity: $K_D = 6.99 \times 10^{-10}$
Isotype: IgG	Full Name: phospholipid transfer protein	$K_{off} = 1.12 \times 10^{-4}$ $K_{on} = 1.60 \times 10^5$
Immunogen Catalog Number: AG0681	Calculated MW: 54.7 kDa, 81 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C .
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00215-1, PLTP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83222-1-PBS. Detection antibody: 83222-2-PBS. Standard: Ag0681. Range: 1.25-20 ng/mL