

For Research Use Only

nectin 4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83223-2-PBS

Basic Information

Catalog Number: 83223-2-PBS	GenBank Accession Number: BC010423	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 81607	CloneNo.: 240009G5
Source: Rabbit	UNIPROT ID: Q96NY8	
Isotype: IgG	Full Name: poliovirus receptor-related 4	
	Calculated MW: 510 aa, 55 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Nectin-4, also known as PVRL4, is a type-I transmembrane glycoprotein that belongs to the nectin subfamily of immunoglobulin-like adhesion molecules that participate in Ca(2+)-independent cell-cell adhesion (PMID: 18703497). The extracellular domain of Nectin-4, which contains three immunoglobulin-like domains (V and two C2-type domains, VCC), can be proteolytically cleaved to release a soluble form. Nectin-4 interacts with afadin via its carboxyl-terminal cytoplasmic sequence and trans-interacts with nectin-1/PRR1 through V domain interaction (PMID: 11544254). It acts as a receptor for measles virus (PMID: 22048310). Nectin-4 is overexpressed in multiple human malignancies and the aberrant expression is correlated with cancer progression and poor prognostic (PMID: 35131876). Nectin-4 can be detected as a soluble form of 43 kDa and a transmembrane form of 66 kDa (PMID: 15784625).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00147-3, nectin 4 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83223-5-PBS. Detection antibody: 83223-2-PBS. Standard: Eg0415. Range: 1.25-20 ng/mL