

For Research Use Only

Parkin Recombinant antibody, PBS Only (Detector)



Catalog Number: 83224-1-PBS

Basic Information

Catalog Number:

83224-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5179

GenBank Accession Number:

BC022014

GeneID (NCBI):

5071

UNIPROT ID:

O60260

Full Name:

Parkinson disease (autosomal recessive, juvenile) 2, parkin

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

240077B8

Applications

Tested Applications:

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00169-1, Human Parkin Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83224-3-PBS. Detection antibody: 83224-1-PBS. Standard: Ag5179. Range: 5-80 ng/mL.



Standard curve of MP00169-2, Human Parkin Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83224-2-PBS. Detection antibody: 83224-1-PBS. Standard: Ag5179. Range: 5-80 ng/mL.