

For Research Use Only

SP8 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83230-3-PBS

Basic Information

Catalog Number: 83230-3-PBS	GenBank Accession Number: BC038669	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 221833	CloneNo.: 240011G4
Source: Rabbit	UNIPROT ID: Q8IXZ3	
Isotype: IgG	Full Name: Sp8 transcription factor	
Immunogen Catalog Number: AG22934	Calculated MW: 508 aa, 51 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00190-1, SP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83230-1-PBS. Detection antibody: 83230-3-PBS. Standard: Ag22934. Range: 0.625-20 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 83230-3-RR against Human SP8 were performed. The affinity constant is 5.44 nM.