

For Research Use Only

PPP1R1A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83231-3-PBS

Basic Information

Catalog Number: 83231-3-PBS	GenBank Accession Number: BC022470	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5502	CloneNo.: 240046G3
Source: Rabbit	UNIPROT ID: Q13522	
Isotype: IgG	Full Name: protein phosphatase 1, regulatory (inhibitor) subunit 1A	
Immunogen Catalog Number: AG3128	Calculated MW: 171 aa, 19 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00197-1, PPP1R1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83231-3-PBS. Detection antibody: 83231-1-PBS. Standard: Ag3128. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00197-2, PPP1R1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83231-3-PBS. Detection antibody: 83231-2-PBS. Standard: Ag3128. Range: 0.625-20 ng/mL.