

For Research Use Only

MTNR1A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83234-3-PBS

Basic Information

Catalog Number: 83234-3-PBS	GenBank Accession Number: BC074947	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4543	CloneNo.: 240005F11
Source: Rabbit	UNIPROT ID: P48039	
Isotype: IgG	Full Name: melatonin receptor 1A	
Immunogen Catalog Number: AG24198	Calculated MW: 350 aa, 39 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00214-2, MTNR1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83234-2-PBS. Detection antibody: 83234-3-PBS. Standard: Ag24198. Range: 1.25-20 ng/mL.