

For Research Use Only

SLC22A10 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83237-5-PBS

Basic Information

Catalog Number:

83237-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34521

GenBank Accession Number:

NM_001039752

GeneID (NCBI):

387775

UNIPROT ID:

Q63ZE4

Full Name:

solute carrier family 22, member 10

Calculated MW:

541 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

240083G6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00221-3, SLC22A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83237-4-PBS. Detection antibody: 83237-5-PBS. Standard: Ag34521. Range: 0.625-20 ng/mL.