

For Research Use Only

USP30 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83240-3-PBS

Basic Information

Catalog Number: 83240-3-PBS	GenBank Accession Number: BC004868	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 84749	CloneNo.: 240111B8
Source: Rabbit	UNIPROT ID: Q70CQ3	
Isotype: IgG	Full Name: ubiquitin specific peptidase 30	
Immunogen Catalog Number: AG7620	Calculated MW: 58.5 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00154-1, USP30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83240-2-PBS. Detection antibody: 83240-3-PBS. Standard: Ag7620. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00154-2, USP30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83240-3-PBS. Detection antibody: 83240-1-PBS. Standard: Ag7620. Range: 0.625-20 ng/mL