

For Research Use Only

Glypican 2 Recombinant antibody, PBS Only (Capture)

Catalog Number:83241-4-PBS



Basic Information

Catalog Number: 83241-4-PBS	GenBank Accession Number: BC027972	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 221914	CloneNo.: 230518F3
Source: Rabbit	UNIPROT ID: Q8N158	Affinity: $K_D < 1.00 \times 10^{-12}$
Isotype: IgG	Full Name: glypican 2	$K_{off} < 1.00 \times 10^{-6}$
Immunogen Catalog Number: AG4041	Calculated MW: 579 aa, 63 kDa	$K_{on} = 7.67 \times 10^4$

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00156-1, GPC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83241-4-PBS. Detection antibody: 83241-3-PBS. Standard: Ag4041. Range: 2.5-80 ng/mL.