

For Research Use Only

# IL3RA/CD123 Recombinant antibody, PBS Only



Catalog Number: 83243-3-PBS

## Basic Information

Catalog Number:

83243-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4452

GenBank Accession Number:

BC035407

GeneID (NCBI):

3563

ENSEMBL Gene ID:

ENSG00000185291

UNIPROT ID:

P26951

Full Name:

interleukin 3 receptor, alpha (low affinity)

Calculated MW:

378 aa, 43 kDa

Purification Method:

Protein A purification

CloneNo.:

240042E5

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HL-60 cells were subjected to SDS PAGE followed by western blot with 83243-3-RR (IL3RA/CD123 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83243-3-PBS in a different storage buffer formulation.