

For Research Use Only

# PRRG2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83244-1-PBS

## Basic Information

Catalog Number: 83244-1-PBS	GenBank Accession Number: BC026032	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5639	CloneNo.: 240074A1
Source: Rabbit	UNIPROT ID: O14669	
Isotype: IgG	Full Name: proline rich Gla (G-carboxyglutamic acid) 2	
Immunogen Catalog Number: AG35095	Calculated MW: 179 aa, 20 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00180-2, PRRG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83244-2-PBS. Detection antibody: 83244-1-PBS. Standard: Ag35095. Range: 2.5-80 ng/mL.