

For Research Use Only

PRRG2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83244-2-PBS

Basic Information

Catalog Number: 83244-2-PBS	GenBank Accession Number: BC026032	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5639	CloneNo.: 240074F12
Source: Rabbit	UNIPROT ID: O14669	
Isotype: IgG	Full Name: proline rich Gla (G-carboxyglutamic acid) 2	
Immunogen Catalog Number: AG35095	Calculated MW: 179 aa, 20 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00180-1, PRRG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83244-2-PBS. Detection antibody: 83244-3-PBS. Standard: Ag35095. Range: 2.5-80 ng/mL.



Cytometric bead array standard curve of MP00180-2, PRRG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83244-2-PBS. Detection antibody: 83244-1-PBS. Standard: Ag35095. Range: 2.5-80 ng/mL.