

For Research Use Only

GJA5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83247-2-PBS

Basic Information

Catalog Number: 83247-2-PBS	GenBank Accession Number: BC013313	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2702	CloneNo.: 240072E10
Source: Rabbit	UNIPROT ID: P36382	
Isotype: IgG	Full Name: gap junction protein, alpha 5, 40kDa	
Immunogen Catalog Number: AG24199	Calculated MW: 40 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00208-1, GJA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83247-3-PBS. Detection antibody: 83247-2-PBS. Standard: Ag24199. Range: 6.25-100 ng/mL