

For Research Use Only

# GJA5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83247-4-PBS

## Basic Information

Catalog Number:

83247-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24199

GenBank Accession Number:

BC013313

GeneID (NCBI):

2702

UNIPROT ID:

P36382

Full Name:

gap junction protein, alpha 5, 40kDa

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

240072F2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00208-2, GJA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83247-4-PBS. Detection antibody: 83247-1-PBS. Standard: Ag24199. Range: 6.25-100 ng/mL