

For Research Use Only

ABHD6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83249-1-PBS

Basic Information

Catalog Number: 83249-1-PBS	GenBank Accession Number: BC001698	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 57406	CloneNo.: 230472C1
Source: Rabbit	UNIPROT ID: Q9BV23	
Isotype: IgG	Full Name: abhydrolase domain containing 6	
Immunogen Catalog Number: AG13968	Calculated MW: 337 aa, 38 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00225-1, ABHD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83249-1-PBS. Detection antibody: 83249-2-PBS. Standard: Ag13968. Range: 3.125-100 ng/mL.