

For Research Use Only

ABHD6 Recombinant antibody, PBS Only



Catalog Number: 83249-3-PBS

Basic Information

Catalog Number: 83249-3-PBS	GenBank Accession Number: BC001698	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57406	CloneNo.: 230472G8
Source: Rabbit	UNIPROT ID: Q9BV23	
Isotype: IgG	Full Name: abhydrolase domain containing 6	
Immunogen Catalog Number: AG13968	Calculated MW: 337 aa, 38 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83249-3-RR (ABHD6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83249-3-PBS in a different storage buffer formulation.