

For Research Use Only

NRG1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83251-1-PBS

Basic Information

Catalog Number: 83251-1-PBS	GenBank Accession Number: BC007675	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3084	CloneNo.: 240067A4
Source: Rabbit	UNIPROT ID: Q02297	
Isotype: IgG	Full Name: neuregulin 1	
Immunogen Catalog Number: AG35114	Calculated MW: 70 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00232-2, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83251-3-PBS. Detection antibody: 83251-1-PBS. Standard: Ag35114. Range: 1.25-20 ng/mL