

For Research Use Only

ND6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83252-2-PBS



Basic Information

Catalog Number:

83252-2-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34496

GenBank Accession Number:

YP_003024037

GeneID (NCBI):

4541

UNIPROT ID:

P03923

Full Name:

NADH dehydrogenase, subunit 6 (complex I)

Calculated MW:

19 kDa

Purification Method:

Protein A purification

CloneNo.:

240146F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00246-2, ND6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83252-4-PBS. Detection antibody: 83252-2-PBS. Standard: Ag34496. Range: 2.5-40 ng/mL.