

For Research Use Only

# TNFRSF11A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83253-2-PBS

## Basic Information

Catalog Number:

83253-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29425

GenBank Accession Number:

NM\_003839

GeneID (NCBI):

8792

UNIPROT ID:

Q9Y6Q6

Full Name:

tumor necrosis factor receptor  
superfamily, member 11a, NFkB  
activator

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

240106E7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00249-1, TNFRSF 11A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83253-2-PBS. Detection antibody: 83253-1-PBS. Standard: Ag29425. Range: 2.5-40 ng/mL.



Cytometric bead array standard curve of MP00249-2, TNFRSF 11A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83253-2-PBS. Detection antibody: 83253-3-PBS. Standard: Ag29425. Range: 2.5-40 ng/mL.