

For Research Use Only

Mouse CD19 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83254-4-PBS

Basic Information

Catalog Number:

83254-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009844.2

GeneID (NCBI):

12478

UNIPROT ID:

P25918

Full Name:

CD19 antigen

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

240135C6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00253-3, Mouse CD19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83254-1-PBS. Detection antibody: 83254-4-PBS. Standard: Eg0907. Range: 5-80 ng/mL.