

For Research Use Only

TMEM40 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83256-3-PBS

Basic Information

| | | |
|--------------------------------------|--|--|
| Catalog Number: 83256-3-PBS | GenBank Accession Number: BC020658 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 55287 | CloneNo.: 240102E6 |
| Source: Rabbit | UNIPROT ID: Q8WWA1 | |
| Isotype: IgG | Full Name: transmembrane protein 40 | |
| Immunogen Catalog Number: AG27882 | | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00163-1, TMEM40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83256-3-PBS. Detection antibody: 83256-1-PBS. Standard: Ag27882. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00163-2, TMEM40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83256-2-PBS. Detection antibody: 83256-3-PBS. Standard: Ag27882. Range: 0.625-20 ng/mL.