

For Research Use Only

SNAP25 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83259-3-PBS

Basic Information

Catalog Number: 83259-3-PBS	GenBank Accession Number: BC010647	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6616	CloneNo.: 240069A3
Source: Rabbit	UNIPROT ID: P60880	
Isotype: IgG	Full Name: synaptosomal-associated protein, 25kDa	
Immunogen Catalog Number: AG6695	Calculated MW: 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00196-1, SNAP25 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83259-3-PBS. Detection antibody: 83259-1-PBS. Standard: Ag6695. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00196-3, SNAP25 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83259-3-PBS. Detection antibody: 83259-2-PBS. Standard: Ag6695. Range: 0.625-20 ng/mL.