

For Research Use Only

PLK1-phospho Recombinant antibody, PBS Only (Detector)



Catalog Number: 83260-3-PBS

Basic Information

Catalog Number:

83260-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005030

GeneID (NCBI):

5347

UNIPROT ID:

P53350

Full Name:

polo-like kinase 1 (Drosophila)

Calculated MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

240092A12

Applications

Tested Applications:

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00238-1, Human PLK1-phospho Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83260-2-PBS. Detection antibody: 83260-3-PBS. Standard: SY00401. Range: 0.156-10 ng/mL