

For Research Use Only

MS4A12 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83261-1-PBS

Basic Information

Catalog Number: 83261-1-PBS	GenBank Accession Number: BC029793	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 54860	CloneNo.: 240105A11
Source: Rabbit	UNIPROT ID: Q9NND0	
Isotype: IgG	Full Name: membrane-spanning 4-domains, subfamily A, member 12	
Immunogen Catalog Number: AG3968	Calculated MW: 267 aa, 28 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00244-1, MS4A12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83261-3-PBS. Detection antibody: 83261-1-PBS. Standard: Ag3968. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00244-2, MS4A12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83261-2-PBS. Detection antibody: 83261-1-PBS. Standard: Ag3968. Range: 0.625-20 ng/mL.