

For Research Use Only

MS4A12 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83261-2-PBS

Basic Information

Catalog Number:

83261-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3968

GenBank Accession Number:

BC029793

GeneID (NCBI):

54860

UNIPROT ID:

Q9N110

Full Name:

membrane-spanning 4-domains,
subfamily A, member 12

Calculated MW:

267 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

240105D5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00244-2, MS4A12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83261-2-PBS. Detection antibody: 83261-1-PBS. Standard: Ag3968. Range: 0.625-20 ng/mL.