

For Research Use Only

# MS4A12 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83261-3-PBS

## Basic Information

Catalog Number:

83261-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3968

GenBank Accession Number:

BC029793

GeneID (NCBI):

54860

UNIPROT ID:

Q9NND0

Full Name:

membrane-spanning 4-domains,  
subfamily A, member 12

Calculated MW:

267 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

240105D10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00244-1, MS4A12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83261-3-PBS. Detection antibody: 83261-1-PBS. Standard: Ag3968. Range: 0.625-20 ng/mL.