

For Research Use Only

EFHD2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83264-1-PBS

Basic Information

Catalog Number: 83264-1-PBS	GenBank Accession Number: BC023611	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 79180	CloneNo.: 240185C3
Source: Rabbit	UNIPROT ID: Q96C19	
Isotype: IgG	Full Name: EF-hand domain family, member D2	
Immunogen Catalog Number: AG35084	Calculated MW: 27 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00269-1, EFHD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83264-1-PBS. Detection antibody: 83264-4-PBS. Standard: Ag35084. Range: 1.25-40 ng/mL



Cytometric bead array standard curve of MP00269-2, EFHD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83264-2-PBS. Detection antibody: 83264-1-PBS. Standard: Ag35084. Range: 1.25-40 ng/mL