

For Research Use Only

DUS2L Recombinant antibody, PBS Only (Capture)

Catalog Number: 83265-1-PBS



Basic Information

Catalog Number: 83265-1-PBS	GenBank Accession Number: BC006527	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 54920	CloneNo.: 240192A6
Source: Rabbit	UNIPROT ID: Q9NX74	
Isotype: IgG	Full Name: dihydrouridine synthase 2-like, SMM1 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG23831		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00261-2, DUS2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83265-3-PBS. Detection antibody: 83265-4-PBS. Standard: Ag23831. Range: 1.25-40 ng/mL