

For Research Use Only

# DUS2L Recombinant antibody, PBS Only (Capture)

Catalog Number: 83265-2-PBS



## Basic Information

<b>Catalog Number:</b> 83265-2-PBS	<b>GenBank Accession Number:</b> BC006527	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 54920	<b>CloneNo.:</b> 240192F8
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9NX74	
<b>Isotype:</b> IgG	<b>Full Name:</b> dihydrouridine synthase 2-like, SMM1 homolog (S. cerevisiae)	
<b>Immunogen Catalog Number:</b> AG23831		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00261-3, DUS2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83265-3-PBS. Detection antibody: 83265-4-PBS. Standard: Ag23831. Range: 1.25-40 ng/mL



Bilayer interferometry (BLI) kinetic assays of 83265-2-RR against Human DUS2L were performed. The affinity constant is 2.66 nM.