

For Research Use Only

BUB3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83266-3-PBS

Basic Information

Catalog Number: 83266-3-PBS	GenBank Accession Number: BC005138	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9184	CloneNo.: 240087D6
Source: Rabbit	UNIPROT ID: O43684	
Isotype: IgG	Full Name: budding uninhibited by benzimidazoles 3 homolog (yeast)	
Immunogen Catalog Number: AG25608	Calculated MW: 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00297-1, BUB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83266-3-PBS. Detection antibody: 83266-4-PBS. Standard: Ag25608. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00297-2, BUB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83266-1-PBS. Detection antibody: 83266-3-PBS. Standard: Ag25608. Range: 1.25-20 ng/mL