

For Research Use Only

CCL17 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83278-2-PBS

Basic Information

Catalog Number: 83278-2-PBS	GenBank Accession Number: BC112066	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 6361	CloneNo.: 230354C8
Source: Rabbit	UNIPROT ID: Q92583	
Isotype: IgG	Full Name: chemokine (C-C motif) ligand 17	
	Calculated MW: 94 aa, 11 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00308-1, CCL17 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83278-1-PBS. Detection antibody: 83278-2-PBS. Standard: Eg0379. Range: 0.625-20 ng/mL