

For Research Use Only

Shh Recombinant antibody

Catalog Number: 83280-2-RR



Basic Information

Catalog Number:

83280-2-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009170.3

GeneID (NCBI):

20423

UNIPROT ID:

Q62226

Full Name:

sonic hedgehog

Calculated MW:

48KD

Observed MW:

19-50 kDa

Purification Method:

Protein A purification

CloneNo.:

240057G4

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : HepG2 cells, C2C12 cells, RAW 264.7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83280-2-RR (SHH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.