

For Research Use Only

MPO Recombinant antibody, PBS Only (Capture)

Catalog Number: 83282-2-PBS



Basic Information

Catalog Number:

83282-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC130476

GeneID (NCBI):

4353

UNIPROT ID:

P05164

Full Name:

myeloperoxidase

Calculated MW:

745 aa, 84 kDa

Purification Method:

Protein A purification

CloneNo.:

240229G2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00321-1, MPO Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83282-2-PBS. Detection antibody: 83282-1-PBS. Standard: Eg0511. Range: 3.13-100 ng/mL