

For Research Use Only

IGFBP2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83283-1-PBS

Basic Information

Catalog Number: 83283-1-PBS	GenBank Accession Number: NM_000597.3	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3485	CloneNo.: 240227F2
Source: Rabbit	UNIPROT ID: P18065	
Isotype: IgG	Full Name: IGF binding protein 2	
	Calculated MW: 35 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00322-1, IGFBP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83283-2-PBS. Detection antibody: 83283-1-PBS. Standard: Eg0823. Range: 0.156-20 ng/mL