

For Research Use Only

# IL6ST Recombinant antibody, PBS Only (Capture)



Catalog Number: 83284-1-PBS

## Basic Information

|                                       |  |   |
|---------------------------------------|--|---|
| <b>Catalog Number:</b><br>83284-1-PBS | <b>GenBank Accession Number:</b><br>NM_002184.3  | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1mg/ml                | <b>GeneID (NCBI):</b><br>3572  | <b>CloneNo.:</b><br>240288A10                         |
| <b>Source:</b><br>Rabbit              | <b>UNIPROT ID:</b><br>P40189-1   |   |
| <b>Isotype:</b><br>IgG                | <b>Full Name:</b><br>interleukin 6 signal transducer<br>(gp130, oncostatin M receptor) |   |
|                                       | <b>Calculated MW:</b><br>104 kDa   |   |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00324-1, gp130 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83284-1-PBS. Detection antibody: 83284-3-PBS. Standard: Eg0654. Range: 0.78-100ng/mL



Cytometric bead array standard curve of MP00324-2, gp130 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83284-1-PBS. Detection antibody: 83284-2-PBS. Standard: Eg0654. Range: 0.78-100ng/mL