

For Research Use Only

# DDX47 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83287-1-PBS

## Basic Information

|                                             |                                                                |                                                       |
|---------------------------------------------|----------------------------------------------------------------|-------------------------------------------------------|
| <b>Catalog Number:</b><br>83287-1-PBS       | <b>GenBank Accession Number:</b><br>BC068009                   | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1mg/ml                      | <b>GeneID (NCBI):</b><br>51202                                 | <b>CloneNo.:</b><br>240193A6                          |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q9H0S4                                   |                                                       |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>DEAD (Asp-Glu-Ala-Asp) box polypeptide 47 |                                                       |
| <b>Immunogen Catalog Number:</b><br>AG11953 | <b>Calculated MW:</b><br>455 aa, 51 kDa                        |                                                       |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00325-1, DDX47 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83287-2-PBS. Detection antibody: 83287-1-PBS. Standard: Ag11953. Range: 0.625-80 ng/mL