

For Research Use Only

DDX47 Recombinant antibody, PBS Only (Capture)



Catalog Number:83287-2-PBS

Basic Information

Catalog Number: 83287-2-PBS	GenBank Accession Number: BC068009	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 51202	CloneNo.: 240193E2
Source: Rabbit	UNIPROT ID: Q9H0S4	
Isotype: IgG	Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 47	
Immunogen Catalog Number: AG11953	Calculated MW: 455 aa, 51 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00325-1, DDX47 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83287-2-PBS. Detection antibody: 83287-1-PBS. Standard: Ag11953. Range: 0.625-80 ng/mL