

For Research Use Only

GAPDHS Recombinant antibody, PBS Only (Detector)



Catalog Number: 83290-2-PBS

Basic Information

Catalog Number:

83290-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4988

GenBank Accession Number:

BC036373

GeneID (NCBI):

26330

UNIPROT ID:

O14556

Full Name:

glyceraldehyde-3-phosphate
dehydrogenase, spermatogenic

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

240277H3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00323-1, GAPDHs Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83290-1-PBS. Detection antibody: 83290-2-PBS. Standard: Ag4988. Range: 3.13-100 ng/mL.