

For Research Use Only

GIN51 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83293-4-PBS

Basic Information

Catalog Number: 83293-4-PBS	GenBank Accession Number: BC012542	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9837	CloneNo.: 240234F10
Source: Rabbit	UNIPROT ID: Q14691	
Isotype: IgG	Full Name: GIN5 complex subunit 1 (Psf1 homolog)	
Immunogen Catalog Number: AG8895	Calculated MW: 196 aa, 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00317-1, GINS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83293-3-PBS. Detection antibody: 83293-4-PBS. Standard: Ag8895. Range: 3.13-100 ng/mL