

For Research Use Only

ECSIT Recombinant antibody, PBS Only (Capture)



Catalog Number: 83295-4-PBS

Basic Information

Catalog Number:

83295-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26898

GenBank Accession Number:

BC008279

GeneID (NCBI):

51295

UNIPROT ID:

Q9BQ95

Full Name:

ECSIT homolog (Drosophila)

Calculated MW:

49 kDa

Purification Method:

Protein A purification

CloneNo.:

240223H8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00315-1, ECSIT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83295-4-PBS. Detection antibody: 83295-1-PBS. Standard: Ag26898. Range: 3.13-100 ng/mL