

For Research Use Only

FGD6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83296-2-PBS

Basic Information

Catalog Number: 83296-2-PBS	GenBank Accession Number: BC013319	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 55785	CloneNo.: 240260C1
Source: Rabbit	UNIPROT ID: Q6ZV73	
Isotype: IgG	Full Name: FYVE, RhoGEF and PH domain containing 6	
Immunogen Catalog Number: AG22211	Calculated MW: 1430 aa, 161 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00314-2, FGD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83296-1-PBS. Detection antibody: 83296-2-PBS. Standard: Ag22211. Range: 3.13-100 ng/mL