

For Research Use Only

FGD6 Recombinant antibody, PBS Only

Catalog Number: 83296-4-PBS



Basic Information

Catalog Number: 83296-4-PBS	GenBank Accession Number: BC013319	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55785	CloneNo.: 240260A6
Source: Rabbit	UNIPROT ID: Q6ZV73	
Isotype: IgG	Full Name: FYVE, RhoGEF and PH domain containing 6	
Immunogen Catalog Number: AG22211	Calculated MW: 1430 aa, 161 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83296-4-RR (FGD6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83296-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83296-4-RR against Human FGD6 were performed. The affinity constant is 10.7 nM.