

For Research Use Only

Phospho-Histone H2A.X (Ser139) Recombinant antibody, PBS Only



Catalog Number: 83307-2-PBS

Basic Information

Catalog Number:

83307-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC013416

GeneID (NCBI):

3014

UNIPROT ID:

P16104

Full Name:

H2A histone family, member X

Calculated MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

5N19

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Non-treated Jurkat cells, and staurosporine treated Jurkat cells were subjected to SDS PAGE followed by western blot with 83307-2-RRR (Phospho-Histone H2A.X (Ser139) antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Alpha Tubulin (66031-1-Ig) antibody as loading control. This data was developed using the same antibody clone with 83307-2-PBS in a different storage buffer formulation.

